

Product Datasheet

Mouse Cytochrome P450 7A1 (CYP7A1) ELISA Kit (orb1206534)

Description	Mouse Cytochrome P450 7A1 (CYP7A1) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q64505
Sensitivity	0.139ng/mL
Entrez	13122
Expiration Date	Please enquire.