

Product Datasheet

Rat Cytochrome P450 2E1 (CYP2E1) ELISA Kit (orb1206535)

Description	Rat Cytochrome P450 2E1 (CYP2E1) ELISA Kit
Reactivity	Rat
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P05182
Sensitivity	1.25ng/mL
Entrez	25086
Expiration Date	Please enquire.