

Product Datasheet

Rat Collagen Type I Alpha 2 (COL1a2) ELISA Kit (orb1206561)

Description	Rat Collagen Type I Alpha 2 (COL1a2) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P02466
Sensitivity	0.55ng/mL
Entrez	84352
Expiration Date	Please enquire.