

Product Datasheet

Rabbit Ciliary Neurotrophic Factor (CNTF) ELISA Kit (orb1206564)

Description	Rabbit Ciliary Neurotrophic Factor (CNTF) ELISA Kit
Reactivity	Rabbit
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P14188
Sensitivity	11.8pg/mL
Entrez	100009117
Expiration Date	Please enquire.