

## Product Datasheet

### Mouse Citrate Synthase (CS) ELISA Kit (orb1206570)

|                          |                                                              |
|--------------------------|--------------------------------------------------------------|
| <b>Description</b>       | Mouse Citrate Synthase (CS) ELISA Kit                        |
| <b>Reactivity</b>        | Mouse                                                        |
| <b>Range</b>             | 1.56-100ng/mL                                                |
| <b>Note</b>              | For research use only                                        |
| <b>Application notes</b> | sandwich                                                     |
| <b>Sample Types</b>      | tissue homogenates, cell lysates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>Q9CZU6</b>                                                |
| <b>Sensitivity</b>       | 0.56ng/mL                                                    |
| <b>Entrez</b>            | <b>12974</b>                                                 |
| <b>Expiration Date</b>   | Please enquire.                                              |