

Product Datasheet

Rat Cholinergic Receptor, Nicotinic, Alpha 1 (CHRNA1) ELISA Kit (orb1206571)

Description	Rat Cholinergic Receptor, Nicotinic, Alpha 1 (CHRNA1) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P25108
Sensitivity	0.122ng/mL
Entrez	79557
Expiration Date	Please enquire.