

Product Datasheet

Mouse Cholinergic Receptor, Muscarinic 3 (CHRM3) ELISA Kit (orb1206572)

Description	Mouse Cholinergic Receptor, Muscarinic 3 (CHRM3) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9ERZ3
Sensitivity	0.105ng/mL
Entrez	12671
Expiration Date	Please enquire.