

Product Datasheet

Mouse Cholinergic Receptor, Muscarinic 2 (CHRM2) ELISA Kit (orb1206573)

Description	Mouse Cholinergic Receptor, Muscarinic 2 (CHRM2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9ERZ4
Sensitivity	0.057ng/mL
Entrez	243764
Expiration Date	Please enquire.