

Product Datasheet

Rabbit Cholesteryl Ester Transfer Protein (CETP) ELISA Kit (orb1206579)

Description	Rabbit Cholesteryl Ester Transfer Protein (CETP) ELISA Kit
Reactivity	Rabbit
Range	93.7-6000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P22687
Sensitivity	42.61pg/mL
Entrez	100327267
Expiration Date	Please enquire.