

Product Datasheet

Mouse Cluster Of Differentiation 83 (CD83) ELISA Kit (orb1206585)

Description	Mouse Cluster Of Differentiation 83 (CD83) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O88324
Sensitivity	35pg/mL
Entrez	12522
Expiration Date	Please enquire.