

Product Datasheet

Mouse Cluster Of Differentiation 4 (CD4) ELISA Kit (orb1206587)

Description	Mouse Cluster Of Differentiation 4 (CD4) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06332
Sensitivity	0.116ng/mL
Entrez	12504
Expiration Date	Please enquire.