

Product Datasheet

Mouse Chemokine C-C-Motif Receptor 7 (CCR7) ELISA Kit (orb1206595)

Description	Mouse Chemokine C-C-Motif Receptor 7 (CCR7) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P47774
Sensitivity	0.054ng/mL
Entrez	12775
Expiration Date	Please enquire.