

Product Datasheet

Mouse Chemokine C-C-Motif Receptor 10 (CCR10) ELISA Kit (orb1206597)

Description	Mouse Chemokine C-C-Motif Receptor 10 (CCR10) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9JL21
Sensitivity	0.051ng/mL
Entrez	12777
Expiration Date	Please enquire.