

Product Datasheet

Rat Calcium/Calmodulin Dependent Protein Kinase II Gamma (CAMK2g) ELISA Kit (orb1206615)

Description	Rat Calcium/Calmodulin Dependent Protein Kinase II Gamma (CAMK2g) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P11730
Sensitivity	0.123ng/mL
Entrez	171140
Expiration Date	Please enquire.