

Product Datasheet

Mouse Breast And Kidney Expressed Chemokine (BRAK) ELISA Kit (orb1206631)

Description	Mouse Breast And Kidney Expressed Chemokine (BRAK) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9WUQ5
Sensitivity	0.054ng/mL
Entrez	57266
Expiration Date	Please enquire.