

Product Datasheet

Rat ATPase, Na⁺/K⁺ Transporting Beta 3 Polypeptide (ATP1b3) ELISA Kit (orb1206647)

Description	Rat ATPase, Na ⁺ /K ⁺ Transporting Beta 3 Polypeptide (ATP1b3) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q63377
Sensitivity	0.071ng/mL
Entrez	25390
Expiration Date	Please enquire.