

## Product Datasheet

# Mouse Ataxia Telangiectasia Mutated (ATM) ELISA Kit (orb1206648)

|                          |                                                                             |
|--------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Ataxia Telangiectasia Mutated (ATM) ELISA Kit                         |
| <b>Reactivity</b>        | Mouse                                                                       |
| <b>Range</b>             | 15.6-1000pg/mL                                                              |
| <b>Note</b>              | For research use only                                                       |
| <b>Application notes</b> | sandwich                                                                    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>Q62388</b>                                                               |
| <b>Sensitivity</b>       | 5.6pg/mL                                                                    |
| <b>Entrez</b>            | <b>11920</b>                                                                |
| <b>Expiration Date</b>   | Please enquire.                                                             |