

Product Datasheet

Mouse Apolipoprotein A5 (APOA5) ELISA Kit (orb1206667)

Description	Mouse Apolipoprotein A5 (APOA5) ELISA Kit
Reactivity	Mouse
Range	3.12-200pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8C7G5
Sensitivity	1.45pg/mL
Entrez	66113
Expiration Date	Please enquire.