

Product Datasheet

Mouse Apolipoprotein A4 (APOA4) ELISA Kit (orb1206668)

Description	Mouse Apolipoprotein A4 (APOA4) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06728
Sensitivity	38pg/mL
Entrez	11808
Expiration Date	Please enquire.