

Product Datasheet

Mouse Angiopoietin Like Protein 4 (ANGPTL4) ELISA Kit (orb1206676)

Description	Mouse Angiopoietin Like Protein 4 (ANGPTL4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Z1P8
Sensitivity	0.065ng/mL
Entrez	57875
Expiration Date	Please enquire.