

Product Datasheet

Mouse Amylase Alpha 2, Pancreatic (AMY2) ELISA Kit (orb1206682)

Description	Mouse Amylase Alpha 2, Pancreatic (AMY2) ELISA Kit
Reactivity	Mouse
Range	9.37-600pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P00688
Sensitivity	4.1pg/mL
Entrez	109959
Expiration Date	Please enquire.