

Product Datasheet

Rat Alanine Aminotransferase (ALT) ELISA Kit (orb1206685)

Description	Rat Alanine Aminotransferase (ALT) ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P25409
Sensitivity	0.28ng/mL
Entrez	81670
Expiration Date	Please enquire.