

Product Datasheet

Mouse Arachidonate-5-Lipoxygenase (ALOX5) ELISA Kit (orb1206687)

Description	Mouse Arachidonate-5-Lipoxygenase (ALOX5) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P48999
Sensitivity	0.058ng/mL
Entrez	11689
Expiration Date	Please enquire.