

Product Datasheet

Goat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit (orb1206724)

Description	Goat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Goat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13296
Sensitivity	5.9pg/mL
Entrez	100861232
Expiration Date	Please enquire.