

Product Datasheet

Equine Interleukin 13 (IL13) ELISA Kit (orb1206740)

Description	Equine Interleukin 13 (IL13) ELISA Kit
Reactivity	Equine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	B6C802
Sensitivity	6.6pg/mL
Entrez	100034113
Expiration Date	Please enquire.