

Product Datasheet

Goat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit (orb1206745)

Description	Goat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Goat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	A6N998
Sensitivity	0.148ng/mL
Entrez	100860789
Expiration Date	Please enquire.