

Product Datasheet

Human Tumor Necrosis Factor Alpha Induced Protein 6 (TNFaIP6) ELISA Kit (orb1206749)

Description	Human Tumor Necrosis Factor Alpha Induced Protein 6 (TNFaIP6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P98066
Sensitivity	0.056ng/mL
Entrez	7130
Expiration Date	Please enquire.