

Product Datasheet

Human 14-3-3 Protein Epsilon (YWHAE) ELISA Kit (orb1206761)

Description	Human 14-3-3 Protein Epsilon (YWHAE) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P62258
Sensitivity	32.9pg/mL
Entrez	7531
Expiration Date	Please enquire.