

Product Datasheet

Human Xeroderma Pigmentosum, Complementation Group G (XPG) ELISA Kit (orb1206765)

Description	Human Xeroderma Pigmentosum, Complementation Group G (XPG) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P28715
Sensitivity	0.042ng/mL
Entrez	2073
Expiration Date	Please enquire.