

Product Datasheet

Human Target Of The Antiproliferative Antibody 1 (TAPA1) ELISA Kit (orb1206804)

Description	Human Target Of The Antiproliferative Antibody 1 (TAPA1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P60033
Sensitivity	0.119ng/mL
Entrez	975
Expiration Date	Please enquire.