

Product Datasheet

Human Superoxide Dismutase 1, Soluble (SOD1) ELISA Kit (orb1206817)

Description	Human Superoxide Dismutase 1, Soluble (SOD1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P00441
Sensitivity	0.063ng/mL
Entrez	6647
Expiration Date	Please enquire.