

Product Datasheet

Human Suppressors Of Cytokine Signaling 3 (SOCS3) ELISA Kit (orb1206818)

Description	Human Suppressors Of Cytokine Signaling 3 (SOCS3) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O14543
Sensitivity	6.4pg/mL
Entrez	9021
Expiration Date	Please enquire.