

Product Datasheet

Human Secretory Leukocyte Peptidase Inhibitor (SLPI) ELISA Kit (orb1206821)

Description	Human Secretory Leukocyte Peptidase Inhibitor (SLPI) ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P03973
Sensitivity	25.8pg/mL
Entrez	6590
Expiration Date	Please enquire.