

Product Datasheet

Human Sialic Acid Binding Ig Like Lectin 8 (SIGLEC8) ELISA Kit (orb1206822)

Description	Human Sialic Acid Binding Ig Like Lectin 8 (SIGLEC8) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	cell lysates.
Uniprot ID	Q9NYZ4
Sensitivity	0.058ng/mL
Entrez	27181
Expiration Date	Please enquire.