

Product Datasheet

Human Sialic Acid Binding Ig Like Lectin 2 (SIGLEC2) ELISA Kit (orb1206823)

Description	Human Sialic Acid Binding Ig Like Lectin 2 (SIGLEC2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P20273
Sensitivity	0.134ng/mL
Entrez	933
Expiration Date	Please enquire.