

Product Datasheet

Human Selenoprotein P1, Plasma (SEPP1) ELISA Kit (orb1206825)

Description	Human Selenoprotein P1, Plasma (SEPP1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P49908
Sensitivity	0.48pg/mL
Entrez	6414
Expiration Date	Please enquire.