

Product Datasheet

Human Phospholipid Scramblase 1 (PLSCR1) ELISA Kit (orb1206867)

Description	Human Phospholipid Scramblase 1 (PLSCR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	015162
Sensitivity	0.067ng/mL
Entrez	5359
Expiration Date	Please enquire.