

Product Datasheet

Human Platelet Factor 4 Variant 1 (PF4V1) ELISA Kit (orb1206880)

Description	Human Platelet Factor 4 Variant 1 (PF4V1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P10720
Sensitivity	0.8ng/mL
Entrez	5197
Expiration Date	Please enquire.