

Product Datasheet

Human Macrophage Scavenger Receptor 1 (MSR1) ELISA Kit (orb1206926)

Description	Human Macrophage Scavenger Receptor 1 (MSR1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P21757
Sensitivity	0.133ng/mL
Entrez	4481
Expiration Date	Please enquire.