

Product Datasheet

Human Membrane Spanning 4 Domains Subfamily A, Member 1 (MS4A1) ELISA Kit (orb1206928)

Description	Human Membrane Spanning 4 Domains Subfamily A, Member 1 (MS4A1) ELISA Kit
Reactivity	Human
Range	1.25-80ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P11836
Sensitivity	0.45ng/mL
Entrez	931
Expiration Date	Please enquire.