

Product Datasheet

Human Matrix Metalloproteinase 28 (MMP28) ELISA Kit (orb1206932)

Description	Human Matrix Metalloproteinase 28 (MMP28) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H239
Sensitivity	0.065ng/mL
Entrez	79148
Expiration Date	Please enquire.