

Product Datasheet

Human Heat Shock 10kDa Protein 1 (HSP10) ELISA Kit (orb1207011)

Description	Human Heat Shock 10kDa Protein 1 (HSP10) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P61604
Sensitivity	0.62ng/mL
Entrez	3336
Expiration Date	Please enquire.