

Product Datasheet

Human High Mobility Group Nucleosome Binding Domain Protein 1 (HMGN1) ELISA Kit (orb1207014)

Description	Human High Mobility Group Nucleosome Binding Domain Protein 1 (HMGN1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P05114
Sensitivity	0.114ng/mL
Entrez	3150
Expiration Date	Please enquire.