

Product Datasheet

Human Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) ELISA Kit (orb1207051)

Description	Human Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P04406
Sensitivity	0.57ng/mL
Entrez	2597
Expiration Date	Please enquire.