

Product Datasheet

Human Growth Associated Protein 43 (GAP43) ELISA Kit (orb1207052)

Description	Human Growth Associated Protein 43 (GAP43) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P17677
Sensitivity	0.062ng/mL
Entrez	2596
Expiration Date	Please enquire.