

Product Datasheet

Human Receptor II For The Fc Region Of Immunoglobulin G (FcγRII) ELISA Kit (orb1207068)

Description	Human Receptor II For The Fc Region Of Immunoglobulin G (FcγRII) ELISA Kit
Reactivity	Human
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P12318
Sensitivity	15.25pg/mL
Entrez	2212
Expiration Date	Please enquire.