

Product Datasheet

Human Elastase 1, Pancreatic (ELA1) ELISA Kit (orb1207095)

Description	Human Elastase 1, Pancreatic (ELA1) ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UNI1
Sensitivity	23.6pg/mL
Entrez	1990
Expiration Date	Please enquire.