

Product Datasheet

Human Defensin Alpha 5, Paneth Cell Specific (DEFa5) ELISA Kit (orb1207105)

Description	Human Defensin Alpha 5, Paneth Cell Specific (DEFa5) ELISA Kit
Reactivity	Human
Range	3.12-200pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q01523
Sensitivity	1.31pg/mL
Entrez	1670
Expiration Date	Please enquire.