

Product Datasheet

Human Cytochrome P450 7A1 (CYP7A1) ELISA Kit (orb1207113)

Description	Human Cytochrome P450 7A1 (CYP7A1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	CYP7A1
Sensitivity	0.34ng/mL
Entrez	1581
Expiration Date	Please enquire.