

Product Datasheet

Human Chemokine C-X-C-Motif Receptor 7 (CXCR7) ELISA Kit (orb1207115)

Description	Human Chemokine C-X-C-Motif Receptor 7 (CXCR7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P25106
Sensitivity	0.062ng/mL
Entrez	57007
Expiration Date	Please enquire.